

gephyrin antibody

Catalog Number: 110943

Product name

gephyrin antibody

IP: 1:200-1:2000

IF: 1:20-1:200


Specificity

Human, Mouse, Rat; other species not tested.

Validations

Antibody description

gephyrin Rabbit Polyclonal antibody. Positive WB detected in human brain tissue, Jurkat cells. Positive IP detected in HEK-293 cells. Positive IF detected in HEK-293 cells. Observed molecular weight by Western-blot: 93 kDa



Preparation

This antibody was obtained by immunization of gephyrin recombinant protein (Accession Number: NM_020806). Purification method: Antigen affinity purified.

human brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:110943(GPHN antibody) at dilution of 1:400

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

ELISA, WB, IP, IF


Dilutions

Recommended Dilution:

WB: 1:500-1:5000



IP Result of anti-GPHN (IP:Catalog No:110943, 4ug; Detection:Catalog No:110943 1:700) with HEK-293 cells lysate 3200ug.



Immunofluorescent analysis of HEK-293 cells using Catalog No:110943(GPHN Antibody) at

gephyrin antibody



Catalog Number: 110943

dilution of 1:50 and Alexa Fluor 488-conjugated

AffiniPure Goat Anti-Rabbit IgG(H+L)