

FOXN2 Polyclonal Antibody



Catalog Number: 163943

Product name

FOXN2 Polyclonal Antibody

Specificity

Human

Antibody description

Polyclonal antibody to FOXN2

Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 172-431 of human FOXN2 (NP_002149.2).

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

IF

Dilutions

IF 1:50 - 1:200

Validations



Immunofluorescence - FOXN2 Polyclonal Antibody

Immunofluorescence analysis of A549 cells using FOXN2 antibody .