

Plasmolipin antibody



Catalog Number: 113943

Product name

Plasmolipin antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

Plasmolipin Rabbit Polyclonal antibody. Positive WB detected in human testis tissue, human kidney tissue. Observed molecular weight by Western-blot: 18kd

Preparation

This antibody was obtained by immunization of Plasmolipin recombinant protein (Accession Number: NM_015993). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:1000-1:10000

Validations



human testis tissue were subjected to SDS PAGE followed by western blot with Catalog No:113943(PLLIP antibody) at dilution of 1:400