

Nofazinlimab (Synonyms: RECOMBINANT PD-1 MAB CS1003)



Catalog Number: 176943

Product name

Nofazinlimab (Synonyms: RECOMBINANT PD-1 MAB CS1003)

Specificity

Human PDCD1 / PD-1 / CD279

Antibody description

Anti-PDCD1 / PD-1 / CD279 Antibody (Nofazinlimab)

Preparation

Recombinant expression and purified from CHO cells.

Formulation

0.1 M Pro-Ac, 20 mM Arg, pH 5.0

Storage

-80°C for 2 years under sterile conditions; -20°C for 1 year under sterile conditions; Avoid repeated freeze-thaw cycles.

Clonality

Monoclonal


Ig Type

IgG4SP

Applications


ELISA

Validations




SDS-PAGE

Anti-PDCD1 / PD-1 / CD279 Antibody (nofazinlimab) on SDS-PAGE under reducing (R) condition. The purity of the protein is greater than 95%.



SEC-HPLC

The purity of Anti-PDCD1 / PD-1 / CD279 Antibody (nofazinlimab) is 98.2%, determined by SEC-HPLC.



ELISA

Immobilized Human PDCD1 / PD-1 / CD279 His Tag, Mouse PDCD1 / PD-1 / CD279 His Tag, Cynomolgus PDCD1 / PD-1 / CD279 His Tag at 2 ug/mL can bind Anti-PDCD1 / PD-1 / CD279 Antibody (nofazinlimab)