

AKR1CL2 antibody



Catalog Number: 107943

Product name

AKR1CL2 antibody

Specificity

human, mouse, rat; other species not tested.

Antibody description

AKR1CL2 Rabbit Polyclonal antibody. Positive WB detected in mouse testis tissue, RAW 264.7 cells, HT-1080 cells. Positive IF detected in HeLa. Observed molecular weight by Western-blot: 36kd

Preparation

This antibody was obtained by immunization of AKR1CL2 recombinant protein (Accession Number: NM_031436). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IF


Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IF: 1:10-1:100

Validations



mouse testis tissue were subjected to SDS PAGE followed by western blot with Catalog No:107943(AKR1CL2 antibody) at dilution of 1:500



Immunofluorescent analysis of HeLa cells, using AKR1CL2 antibody Catalog No:107943 at 1:25 dilution and Rhodamine-labeled goat anti-rabbit IgG (red). Blue pseudocolor = DAPI (fluorescent DNA dye).