

Anti-DACH2 antibody



Catalog Number: 141243

Product name

Anti-DACH2 antibody

Specificity

Human DACH2

Antibody description

Rabbit Polyclonal to Human DACH2

Preparation

Produced in rabbits immunized with a synthetic peptide corresponding to the center region of the Human DACH2, and purified by antigen affinity chromatography.

Formulation

0.2 µm filtered solution in PBS

Storage

This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C.

Preservative-Free.

Sodium azide is recommended to avoid contamination (final concentration 0.05%-0.1%). It

is toxic to cells and should be disposed of properly. Avoid repeated freeze-thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

WB

Dilutions

WB: 5-20 µg/ml

Validations



Human DACH2 Western blot (WB) 14928