

PARD6A Polyclonal Antibody



Catalog Number: 161243

Product name

PARD6A Polyclonal Antibody

Specificity

Rat

Antibody description

Polyclonal antibody to PARD6A

Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 256-345 of human PARD6A (NP_001032358.1).

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

WB

Dilutions

WB 1:500 - 1:2000

Validations



Western blot - PARD6A Polyclonal Antibody

Western blot analysis of extracts of rat skeletal muscle, using PARD6A antibody at 1:3000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 15s.