

# Human EDNRB / Endothelin B Receptor (Fc Tag) recombinant protein



Catalog Number: 502243

## General Information

### Gene Name Synonym

Endothelin receptor non-selective type

### Protein Construction

A DNA sequence encoding the human EDNRB (NP\_000106.1) (Met1-Lys101) was expressed with the Fc region of mouse IgG1 at the C-terminus.

### Organism

Human

### Expression Host

Human Cells

## QC Testing

### Purity

> 85 % as determined by SDS-PAGE

### Endotoxin

< 1.0 EU per µg of the protein as determined by the LAL method

### Stability

Samples are stable for up to twelve months from date of receipt at -70°C

### Predicted N terminal

Glu 27

### Molecular Mass

The recombinant human EDNRB/mFc comprises

309 amino acids and has a predicted molecular mass of 34.4 kDa. The apparent molecular mass of the protein is approximately 47.5 and 43.7 kDa in SDS-PAGE under reducing conditions.

### Formulation

Lyophilized from sterile PBS, pH 7.4.

1. 5 % trehalose and mannitol are added as protectants before lyophilization.

2. Please contact us for any concerns or special requirements.

## Usage Guide


### Storage

Store it under sterile conditions at -20°C to -80°C. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

### Reconstitution

Adding sterile water, prepare a stock solution of 0.25 mg/ml. Concentration is measured by UV-Vis.

### SDS-PAGE



Human EDNRB / Endothelin B Receptor Protein (Fc Tag) SDS-PAGE