

Anti-CDK5R2 antibody



Catalog Number: 176243

Product name

Anti-CDK5R2 antibody

Specificity

Human

Antibody description

Rabbit polyclonal antibody to CDK5R2

Preparation

This antigen of this antibody was the antiserum was produced against synthesized peptide derived from human cdk5r2. aa range:81-130

Formulation

Liquid, Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Storage

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB, IHC-P

Dilutions

WB: 1:500-1:2000

IHC-P: 1:100-1:300

Validations




Fig1.; Western Blot analysis of various cells using p39 Polyclonal Antibody cells nucleus extracted by Minute TM Cytoplasmic and Nuclear Fractionation kit (SC-003, Inventibiotec, MN, USA).