

# CEACAM21 antibody



Catalog Number: 109243

## Product name

CEACAM21 antibody

## Specificity

Human, Mouse; other species not tested.

## Antibody description

CEACAM21 Rabbit Polyclonal antibody. Positive WB detected in HeLa cells, mouse liver tissue. Positive IF detected in Hela cells. Observed molecular weight by Western-blot: 32,22kd

## Preparation

This antibody was obtained by immunization of CEACAM21 recombinant protein (Accession Number: XM\_047439610). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IF


## Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IF: 1:10-1:100

## Validations



HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:109243(CEACAM21 antibody) at dilution of



1:500 Immunofluorescent analysis of Hela cells, using CEACAM21 antibody Catalog No:109243 at 1:25 dilution and Rhodamine-labeled goat anti-rabbit IgG (red).