

GRM5 Monoclonal Antibody



Catalog Number: 160643

Product name

GRM5 Monoclonal Antibody

Specificity

Human, Mouse, Rat

Antibody description

Monoclonal antibody to GRM5

Preparation

Antigen: Recombinant protein of human GRM5

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Monoclonal

Ig Type

Rabbit IgG

Applications

WB, IHC, IF


Dilutions

WB 1:500 - 1:2000

IHC 1:50 - 1:200


IF 1:50 - 1:200

Validations




Immunohistochemistry - GRM5 Monoclonal Antibody

Immunohistochemistry of paraffin-embedded mouse brain using GRM5 antibody at dilution of 1:100 (40x lens).



Immunohistochemistry - GRM5 Monoclonal Antibody

Immunohistochemistry of paraffin-embedded rat brain using GRM5 antibody at dilution of 1:100




(40x lens).

Immunofluorescence - GRM5 Monoclonal Antibody

GRM5 Monoclonal Antibody

Catalog Number: 160643


Immunofluorescence analysis of SHG-44 cells using GRM5 antibody .



Immunofluorescence - GRM5 Monoclonal Antibody

Immunofluorescence analysis of JAR cells using

GRM5 antibody .



Immunofluorescence - GRM5 Monoclonal Antibody

Immunofluorescence analysis of PC-12 cells using GRM5 antibody .