

RNF125 Polyclonal Antibody



Catalog Number: 165643

Product name

RNF125 Polyclonal Antibody

Specificity

Human

Antibody description

Polyclonal antibody to RNF125

Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 1-232 of human RNF125 (NP_060301.2).

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

IF

Dilutions

IF 1:50 - 1:100

Validations



Immunofluorescence - RNF125 Polyclonal Antibody

Immunofluorescence analysis of A549 cells using RNF125 antibody . Blue: DAPI for nuclear staining.