

SPRR1B antibody

Catalog Number: 115643

Product name

SPRR1B antibody

Specificity

Human; other species not tested.

Antibody description

SPRR1B Rabbit Polyclonal antibody. Positive IHC detected in human skin tissue.

Preparation

This antibody was obtained by immunization of SPRR1B recombinant protein (Accession Number: NM_003125). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, IHC

Dilutions


Recommended Dilution:

IHC: 1:20-1:200

Validations



Immunohistochemistry of paraffin-embedded human skin tissue slide using Catalog No:115643 (SPRR1B Antibody) at dilution of 1:50 (under 40x lens)



Immunohistochemistry of paraffin-embedded human skin tissue slide using Catalog No:115643 (SPRR1B Antibody) at dilution of 1:50 (under 40x lens)