

TRIM24 antibody



Catalog Number: 107643

Product name

TRIM24 antibody

Immunogen

[Human TRIM24 Recombinant protein \(His tag\)](#)

Specificity

Human, Mouse; other species not tested.

Antibody description

TRIM24 Mouse Monoclonal antibody. Positive WB detected in NIH/3T3 cells, HEK-293 cells, HeLa cells. Observed molecular weight by Western-blot: 116 kDa

Preparation

This antibody was obtained by immunization of TRIM24 recombinant protein (Accession Number: NM_015905). Purification method: Protein G purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Monoclonal

Ig Type

Mouse IgG1


Applications

ELISA, WB

Dilutions

Recommended Dilution:
WB: 1:500-1:5000

Validations



NIH/3T3 cells were subjected to SDS PAGE followed by western blot with Catalog No:107643(TRIM24 Antibody) at dilution of 1:1000