

Anti-BLNK Rabbit antibody



Catalog Number: 168643

Product name

Anti-BLNK Rabbit antibody

Specificity

Human, Mouse, Rat

Antibody description

BLNK Rabbit monoclonal antibody

Preparation

Antigen: A synthetic peptide of human BLNK

Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

Storage

Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Monoclonal

Ig Type

Rabbit IgG

Applications


WB, IP

Dilutions

WB: 1/1000

IP: 1/20

Validations



Western blot detection of BLNK in K562, C6, 3T3, HeLa cell lysates using BLNK Rabbit mAb(1:1000 diluted). Predicted band size: 51kDa. Observed band size: 70kDa.