

# Oligophrenin 1 antibody



Catalog Number: 113343

## Product name

Oligophrenin 1 antibody

WB: 1:500-1:5000

IP: 1:200-1:2000

## Specificity


Human, Mouse; other species not tested.

IHC: 1:20-1:200

## Antibody description

Oligophrenin 1 Rabbit Polyclonal antibody. Positive IHC detected in human colon tissue, human stomach cancer tissue. Positive WB detected in mouse brain tissue, mouse skeletal muscle tissue. Positive IP detected in mouse brain tissue. Observed molecular weight by Western-blot: 91 kDa

## Validations




## Preparation

This antibody was obtained by immunization of Oligophrenin 1 recombinant protein (Accession Number: BC059393). Purification method: Antigen affinity purified.

mouse brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:113343(OPHN1 antibody) at dilution of 1:500

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.



## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG


## Applications

ELISA, WB, IHC, IP

## Dilutions

Recommended Dilution:

IP Result of anti-OPHN1 (IP:Catalog No:113343, 4ug; Detection:Catalog No:113343 1:500) with mouse brain tissue lysate 5200ug.




# Oligophrenin 1 antibody

Catalog Number: 113343

---

Immunohistochemistry of paraffin-embedded human colon tissue slide using Catalog No:113343(OPHN1 Antibody) at dilution of 1:50 (under 10x lens)



Immunohistochemistry of paraffin-embedded human colon tissue slide using Catalog No:113343(OPHN1 Antibody) at dilution of 1:50 (under 40x lens)