

# IL37 antibody

Catalog Number: 111743

## Product name

IL37 antibody

## Specificity

Human, Mouse; other species not tested.

## Antibody description

IL37 Rabbit Polyclonal antibody. Positive WB detected in Transfected HEK-293 cells, mouse lung tissue.

## Preparation

This antibody was obtained by immunization of IL37 recombinant protein (Accession Number: XM\_011510963). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:200-1:2000

## Validations



Transfected HEK-293 cells were subjected to SDS PAGE followed by western blot with Catalog No:111743(IL1F7 Antibody) at dilution of 1:500