

IL37 antibody



Catalog Number: 111743

Product name

Polyclonal

IL37 antibody

Ig Type

Specificity

Rabbit IgG

Human, Mouse; other species not tested.

Applications

Antibody description

ELISA, WB

IL37 Rabbit Polyclonal antibody. Positive WB detected in Transfected HEK-293 cells, mouse lung tissue.

Dilutions

Recommended Dilution:

Preparation


WB: 1:200-1:2000

This antibody was obtained by immunization of IL37 recombinant protein (Accession Number: XM_011510963). Purification method: Antigen affinity purified.

Validations

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.



Storage

Store at -20°C. DO NOT ALIQUOT

Transfected HEK-293 cells were subjected to SDS PAGE followed by western blot with Catalog No:111743(IL1F7 Antibody) at dilution of 1:500

Clonality