

RNF151 antibody

Catalog Number: 114743

Product name

RNF151 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

RNF151 Rabbit Polyclonal antibody. Positive WB detected in rat testis tissue, human testis tissue, mouse testis tissue. Positive IHC detected in human testis tissue. Observed molecular weight by Western-blot: 26 kDa

Preparation

This antibody was obtained by immunization of RNF151 recombinant protein (Accession Number: NM_174903). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC


Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IHC: 1:20-1:200


Validations



rat testis tissue were subjected to SDS PAGE followed by western blot with Catalog No:114743(RNF151 Antibody) at dilution of 1:600



Immunohistochemistry of paraffin-embedded human testis tissue slide using Catalog No:114743(RNF151 Antibody) at dilution of 1:50 (under 10x lens)



Immunohistochemistry of paraffin-embedded human testis tissue slide using Catalog No:114743(RNF151 Antibody) at dilution of 1:50

RNF151 antibody



Catalog Number: 114743

(under 40x lens)