

Anti-RAP1A/RAP1B Rabbit antibody



Catalog Number: 166743

Product name

Anti-RAP1A/RAP1B Rabbit antibody

Specificity

Human, Mouse, Rat

Antibody description

RAP1A/RAP1B Rabbit polyclonal antibody

Preparation

Antigen: A synthetic peptide of human RAP1A/B

Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

Storage

Store at -20°C. Stable for 12 months from date of receipt.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB, FC, IP


Dilutions

WB: 1/2000-1/10000

FC: 1/20

IP: 1/20

Validations



Western blot detection of RAP1A/B in K562, Rat Brain, C6, 3T3, HeLa cell lysates using RAP1A/B Rabbit pAb(1:1000 diluted). Predicted band size: 21kDa. Observed band size: 21kDa.