

EXO1 Polyclonal Antibody



Catalog Number: 165053

Product name

EXO1 Polyclonal Antibody

Specificity

Human

Antibody description

Polyclonal antibody to EXO1

Preparation

Antigen: A synthetic peptide of human EXO1

Formulation

PBS with 0.02% sodium azide, pH7.3.

Storage

Store at 4°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


WB, IF

Dilutions

WB 1:500 - 1:1000

IF 1:20 - 1:50

Validations



Western blot - EXO1 Polyclonal Antibody

Western blot analysis of extracts of NCI-H460 cells, using EXO1 antibody. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST.