

# ZHX2 antibody

Catalog Number: 117053

## Product name

ZHX2 antibody

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

ZHX2 Rabbit Polyclonal antibody. Positive IP detected in mouse brain tissue. Positive WB detected in human heart tissue, HeLa cells, Jurkat cells, mouse brain tissue, mouse lung tissue, mouse ovary tissue, mouse skeletal muscle tissue. Observed molecular weight by Western-blot: 92-100 kDa

## Preparation

This antibody was obtained by immunization of ZHX2 recombinant protein (Accession Number: XM\_011516931). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IP


## Dilutions

Recommended Dilution:


WB: 1:1000-1:10000

IP: 1:200-1:2000

## Validations



human heart tissue were subjected to SDS PAGE followed by western blot with Catalog No:117053(ZHX2 antibody) at dilution of 1:800



IP Result of anti-ZHX2 (IP:Catalog No:117053, 4ug; Detection:Catalog No:117053 1:800) with mouse brain tissue lysate 2800ug.