

Anti-C3 Rabbit antibody

Catalog Number: 169053

Product name

Anti-C3 Rabbit antibody

Specificity

Human

Antibody description

C3 Rabbit monoclonal antibody

Preparation

Antigen: A synthetic peptide of human C3

Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

Storage

Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Monoclonal

Ig Type

Rabbit IgG


Applications

WB

Dilutions

WB: 1/1000

Validations



Western blot detection of C3 in K562, HeLa cell lysates using C3 Rabbit mAb(1:1000 diluted). Predicted band size: 187kDa. Observed band size: 120kDa.