

FAH antibody

Catalog Number: 110453

Product name

FAH antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

FAH Rabbit Polyclonal antibody. Positive WB detected in mouse liver tissue, human liver tissue, rat liver tissue. Observed molecular weight by Western-blot: 40-45kd

Preparation

This antibody was obtained by immunization of FAH recombinant protein (Accession Number: NM_000137). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:2000

Validations



mouse liver tissue were subjected to SDS PAGE followed by western blot with Catalog No:110453(FAH antibody) at dilution of 1:500