

# Human TNFRSF17/BCMA?CD269? (Fc Tag) recombinant protein



Catalog Number: 502853

## General Information

### Gene Name Synonym

B-cell maturation protein

### Protein Construction

A DNA sequence encoding the human TNFRSF17 (NP\_001183.2) (Met1-Ala54) was expressed with the Fc region of rabbit IgG at the N-terminus.

### Organism

Human

### Expression Host

Human Cells

## QC Testing

### Purity

> 95 % as determined by SDS-PAGE.

### Endotoxin

< 1.0 EU per  $\mu\text{g}$  protein as determined by the LAL method.

### Stability

Samples are stable for up to twelve months from date of receipt at  $-70^{\circ}\text{C}$

### Predicted N terminal

Ser

### Molecular Mass

The recombinant human TNFRSF17 consists 303 amino acids and predicts a molecular mass of 33.3 kDa.

### Formulation

Lyophilized from sterile PBS, pH 7.4.

1. 5 % trehalose and mannitol are added as protectants before lyophilization.

2. Please contact us for any concerns or special requirements.

## Usage Guide

### Storage


Store it under sterile conditions at  $-20^{\circ}\text{C}$  to  $-80^{\circ}\text{C}$ .

It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

### Reconstitution

Adding sterile water, prepare a stock solution of 0.25 mg/ml. Concentration is measured by UV-Vis.

### SDS-PAGE



Human TNFRSF17 / BCMA / CD269 Protein (Fc Tag)