

PKC zeta antibody



Catalog Number: 113853

Product name

PKC zeta antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

PKC zeta Rabbit Polyclonal antibody. Positive WB detected in HEK-293 cells, NIH/3T3 cells. Positive IHC detected in human breast cancer tissue. Observed molecular weight by Western-blot: 68 kDa

Preparation

This antibody was obtained by immunization of PKC zeta recombinant protein (Accession Number: NM_002744). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, IHC, WB


Dilutions

Recommended Dilution:


WB: 1:200-1:1000

IHC: 1:20-1:200

Validations



Immunohistochemical of paraffin-embedded human breast cancer using Catalog No:113853(PRKCZ antibody) at dilution of 1:100



(under 10x lens)

HEK-293 cells were subjected to SDS PAGE followed by western blot with Catalog No:113853(PRKCZ Antibody) at dilution of 1:300