

WDFY4 antibody



Catalog Number: 116853

Product name

WDFY4 antibody

Specificity

human, mouse, rat; other species not tested.

Antibody description

WDFY4 Rabbit Polyclonal antibody. Positive WB detected in HeLa cells. Observed molecular weight by Western-blot: 73kd

Preparation

This antibody was obtained by immunization of WDFY4 recombinant protein (Accession Number: BC032420). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:500-1:5000

Validations



HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:116853(WDFY4 antibody) at dilution of 1:500