

General Information

Protein Construction

A DNA sequence encoding the extracellular domain of mouse MARCO (NP_034896.1) (Gln 70-Ser 518) was expressed, fused with a polyhistidine tag at the N-terminus.

Organism

Mouse

Expression Host

Human Cells

QC Testing

Purity

> 92 % as determined by SDS-PAGE

Endotoxin

< 1.0 EU per μg of the protein as determined by the LAL method

Stability

Samples are stable for up to twelve months from date of receipt at -70°C

Predicted N terminal

His

Molecular Mass

The secreted recombinant mouse MARCO consists of 465 amino acids and has a calculated molecular mass of 47.3 kDa. As a result of glycosylation, the

recombinant protein migrates as an approximately 55-60 kDa band in SDS-PAGE under reducing conditions.

Formulation

Lyophilized from sterile PBS, pH 7.4

1. 5 % trehalose and mannitol are added as protectants before lyophilization.

2. Please contact us for any concerns or special requirements.

Usage Guide


Storage

Store it under sterile conditions at -20°C to -80°C . It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

Reconstitution

Adding sterile water, prepare a stock solution of 0.25 mg/ml. Concentration is measured by UV-Vis.

SDS-PAGE



Mouse MARCO Protein (His Tag) SDS-PAGE