

NFIX-Specific antibody



Catalog Number: 113153

Product name

NFIX-Specific antibody

Specificity

human, mouse, rat; other species not tested.

Antibody description

NFIX-Specific Rabbit Polyclonal antibody. Positive WB detected in mouse brain tissue, HeLa cells. Observed molecular weight by Western-blot: 23kd

Preparation

This antibody was obtained by immunization of Peptide (Accession Number: NM_002501). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:2000

Validations



mouse brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:113153(NFIX-Specific antibody) at dilution of 1:500