

GATA4 Monoclonal Antibody



Catalog Number: 164153

Product name

GATA4 Monoclonal Antibody

Specificity

Human, Mouse, Rat

Antibody description

Monoclonal antibody to GATA4

Preparation

Antigen: A synthetic peptide of human GATA4

Formulation

PBS with 0.02% sodium azide, pH7.3.

Storage

Store at 4°C. Avoid freeze / thaw cycles.

Clonality

Monoclonal

Ig Type

Mouse IgG

Applications

WB, IHC, IF, IP

Dilutions


WB 1:200 - 1:500

IHC 1:50 - 1:200

IF 1:50 - 1:200


IP 1:50 - 1:100

Validations



Western blot - GATA4 Monoclonal Antibody

Western blot analysis of extracts of mouse liver,



using GATA4 antibody .

Immunofluorescence - GATA4 Monoclonal Antibody

Immunofluorescence analysis of HeLa cells using GATA4 antibody . Blue: DAPI for nuclear staining.