

# HLA-DQA2 Polyclonal Antibody



Catalog Number: 166153

## Product name

HLA-DQA2 Polyclonal Antibody

## Specificity

Mouse, Rat

## Antibody description

Polyclonal antibody to HLA-DQA2

## Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 24-125 of human HLA-DQA2 (NP\_064440.1).

## Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

## Storage

Store at -20°C. Avoid freeze / thaw cycles.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG


## Applications

WB

## Dilutions

WB 1:500 - 1:1000

## Validations



Western blot - HLA-DQA2 Polyclonal Antibody

Western blot analysis of extracts of various cell lines, using HLA-DQA2 antibody at 1:3000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 90s.