

Anti-OR10T2 antibody



Catalog Number: 176153

Product name

Anti-OR10T2 antibody

Specificity

Human

Antibody description

Rabbit polyclonal antibody to OR10T2

Preparation

This antigen of this antibody was the antiserum was produced against synthesized peptide derived from human or10t2. aa range:231-280

Formulation

Liquid, Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Storage

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB

Dilutions

WB: 1:500-1:2000

Validations




Fig1.; Western Blot analysis of 293T-UV cells using Olfactory receptor 10T2 Polyclonal Antibody diluted at 1: 500