

CHD4 antibody



Catalog Number: 107153

Product name

CHD4 antibody

Specificity

Human; other species not tested.

Antibody description

CHD4 Mouse Monoclonal antibody. Positive IF detected in HEK-293 cells. Positive WB detected in HEK-293 cells. Observed molecular weight by Western-blot: 260 kDa

Preparation

This antibody was obtained by immunization of CHD4 recombinant protein (Accession Number: XM_006718958). Purification method: Protein G purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Monoclonal

Ig Type

Mouse IgG1

Applications

ELISA, WB, IF


Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IF: 1:50-1:500

Validations



HEK-293 cells were subjected to SDS PAGE followed by western blot with (CHD4 Antibody) at dilution of 1:500



Immunofluorescent analysis of (10% Formaldehyde) fixed HEK-293 cells using Catalog No:107153(CHD4 Antibody) at dilution of 1:100 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Mouse IgG(H+L)