

# ZNF271 antibody



Catalog Number: 117153

## Product name

ZNF271 antibody

## Specificity

Human, Mouse; other species not tested.

## Antibody description

ZNF271 Rabbit Polyclonal antibody. Positive WB detected in mouse pancreas tissue. Observed molecular weight by Western-blot: 48 kDa

## Preparation

This antibody was obtained by immunization of ZNF271 recombinant protein (Accession Number: BC017710). Purification method: Antigen Affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:500-1:5000

## Validations



mouse pancreas tissue were subjected to SDS PAGE followed by western blot with Catalog No:117153(ZNF271 antibody) at dilution of 1:1000