

# Anti-Phospho-Tyrosine (6G6) Mouse antibody



Catalog Number: 168153

## Product name

Anti-Phospho-Tyrosine (6G6) Mouse antibody

## Specificity

ALL

## Antibody description

Phospho-Tyrosine (6G6) Mouse monoclonal antibody

## Preparation

Antigen: Purified Protein

## Formulation

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

## Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

## Clonality

Monoclonal

## Ig Type

Mouse IgG1

## Applications

WB, IHC-P


## Dilutions

WB: 1/500 - 1/2000


IHC: 1/100 - 1/300

ELISA: 1/40000

## Validations



Western blot analysis of 1) Jurkat, 2) 293T, 3) Rat Brain Tissue, 4) Mouse Brain Tissue with Phosphotyrosine Mouse mAb diluted at 1:2,000



Immunohistochemical analysis of paraffin-embedded Human Breast Carcinoma using Phosphotyrosine Mouse mAb diluted at 1:200.