

# GFRA2 antibody

Catalog Number: 110953

## Product name

GFRA2 antibody

## Specificity

Human, Mouse; other species not tested.

## Antibody description

GFRA2 Rabbit Polyclonal antibody. Positive WB detected in human placenta tissue. Positive IHC detected in human brain tissue, human placenta tissue. Observed molecular weight by Western-blot: 41-45kd

## Preparation

This antibody was obtained by immunization of Peptide (Accession Number: NM\_001495). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IHC


## Dilutions

Recommended Dilution:


WB: 1:200-1:1000

IHC: 1:20-1:200


## Validations



human placenta tissue were subjected to SDS PAGE followed by western blot with Catalog No:110953(GFRA2 antibody) at dilution of 1:200



Immunohistochemistry of paraffin-embedded human brain tissue slide using Catalog No:110953(GFRA2 Antibody) at dilution of 1:50 (under 10x lens)



Immunohistochemistry of paraffin-embedded human brain tissue slide using Catalog No:110953(GFRA2 Antibody) at dilution of 1:50 (under 40x lens)