

TET2 antibody



Catalog Number: 115953

Product name

TET2 antibody

Immunogen

[Human TET2 Recombinant protein \(GST tag\)](#)

Specificity

Human, Mouse; other species not tested.

Antibody description

TET2 Rabbit Polyclonal antibody. Positive IP detected in mouse brain tissue. Positive WB detected in HL-60 cells, A549 cells, HEK-293 cells, human brain tissue. Observed molecular weight by Western-blot: 130kd,133kd,224kd

Preparation

This antibody was obtained by immunization of TET2 recombinant protein (Accession Number: BC150180). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, IP, WB


Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IP: 1:200-1:2000

Validations



HL-60 cells were subjected to SDS PAGE followed by western blot with Catalog No:115953(TET2 antibody) at dilution of 1:500



IP Result of anti-TET2 (IP:Catalog No:115953, 3ug; Detection:Catalog No:115953 1:500) with mouse brain tissue lysate 4000ug.