

# Begelomab (Synonyms: BEGEDINA)



Catalog Number: 176953

## Product name

Begelomab (Synonyms: BEGEDINA)

## Specificity

Human DPP4 / CD26

## Antibody description

Anti-DPP4 / CD26 Antibody (Begelomab)

## Preparation

Recombinant expression and purified from CHO cells.

## Formulation

0.1 M Pro-Ac, 20 mM Arg, pH 5.0

## Storage

-80°C for 2 years under sterile conditions; -20°C for 1 year under sterile conditions; Avoid repeated freeze-thaw cycles.

## Clonality

Monoclonal


## Ig Type

IgG2SA

## Applications


ELISA

## Validations



## SDS-PAGE

Anti-DPP4 / CD26 Antibody (begelomab) on SDS-PAGE under reducing (R) condition. The purity of the protein is greater than 95%.



## SEC-HPLC

The purity of Anti-DPP4 / CD26 Antibody (begelomab) is 98.65%, determined by SEC-HPLC.