

ZNF174 antibody

Catalog Number: 116953

Product name

ZNF174 antibody

IHC: 1:20-1:200

IF: 1:10-1:100

Specificity

human; other species not tested.

Antibody description

ZNF174 Rabbit Polyclonal antibody. Positive WB detected in BxPC-3 cells, HeLa cells, HepG2 cells. Positive IHC detected in human testis. Positive IF detected in Hela. Observed molecular weight by Western-blot: 46kd

Preparation

This antibody was obtained by immunization of ZNF174 recombinant protein (Accession Number: NM_001032292). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB, IHC, IF

Dilutions


Recommended Dilution:


WB: 1:500-1:5000

Validations



BxPC-3 cells were subjected to SDS PAGE followed by western blot with Catalog No:116953(ZNF174 antibody) at dilution of 1:500




Immunohistochemical of paraffin-embedded human testis using Catalog No:116953(ZNF174 antibody) at dilution of 1:100 (under 10x lens) 

Immunohistochemical of paraffin-embedded human testis using Catalog No:116953(ZNF174 antibody) at dilution of 1:100 (under 40x lens)

ZNF174 antibody

Catalog Number: 116953



Immunofluorescent analysis of HeLa cells, using ZNF174 antibody Catalog No:116953 at 1:25 dilution and Rhodamine-labeled goat anti-rabbit IgG (red).