

# Anti-Cdc34 Rabbit antibody



Catalog Number: 167953

## Product name

Anti-Cdc34 Rabbit antibody

## Specificity

Human, Mouse, Rat

## Antibody description

Cdc34 Rabbit monoclonal antibody

## Preparation

Antigen: Recombinant protein of human Cdc34

## Formulation

50nM Tris-Glycine(pH 7.4), 0.15M NaCl,  
40%Glycerol, 0.01% sodium azide and 0.05% BSA

## Storage

Store at 4°C short term. Store at -20°C long term.  
Avoid freeze / thaw cycle.

## Clonality

Monoclonal

## Ig Type

Rabbit IgG

## Applications

WB, ICC/IF, FC, IP

## Dilutions


WB: 1/1000

ICC/IF: 1/20

FC: 1/20

IP: 1/20

## Validations



Western blot detection of Cdc34 in HL-60, U251, U87-MG cell lysates using Cdc34 Rabbit mAb(1:1000 diluted). Predicted band size:27kDa. Observed band size:32kDa.