

BAF60a Polyclonal Antibody



Catalog Number: 166253

Product name

BAF60a Polyclonal Antibody

Specificity

Human, Mouse, Rat

Antibody description

Polyclonal antibody to BAF60a

Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 1-120 of mouse BAF60a (NP_114030.2).

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


WB, IHC

Dilutions

WB 1:500 - 1:2000


IHC 1:50 - 1:200

Validations




Western blot - BAF60a Polyclonal Antibody

Western blot analysis of extracts of various cell lines, using Smarcd1 antibody at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 15s.



Immunohistochemistry - BAF60a Polyclonal Antibody

Immunohistochemistry of paraffin-embedded rat brain using Smarcd1 antibody (40x lens).




BAF60a Polyclonal Antibody



Catalog Number: 166253

Immunohistochemistry - BAF60a Polyclonal Antibody


Immunohistochemistry of paraffin-embedded human lung cancer using Smarcd1 antibody (40x



lens).

Immunohistochemistry - BAF60a Polyclonal Antibody


Immunohistochemistry of paraffin-embedded mouse brain using Smarcd1 antibody (40x lens).



Immunohistochemistry - BAF60a Polyclonal Antibody

Immunohistochemistry of paraffin-embedded rat


brain using Smarcd1 antibody at dilution of 1:100



(40x lens).

Immunohistochemistry - BAF60a Polyclonal Antibody

Immunohistochemistry of paraffin-embedded human lung cancer using Smarcd1 antibody at dilution of 1:100 (40x lens).



Immunohistochemistry - BAF60a Polyclonal Antibody

Immunohistochemistry of paraffin-embedded mouse brain using Smarcd1 antibody at dilution of 1:100 (40x lens).