

Anti-Cdk4 Rabbit antibody



Catalog Number: 167253

Product name

Anti-Cdk4 Rabbit antibody

Specificity

Human

Antibody description

Cdk4 Rabbit polyclonal antibody

Preparation

Antigen: A synthetic peptide of human Cdk4

Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

Storage

Store at -20°C. Stable for 12 months from date of receipt.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB, ICC/IF, FC


Dilutions

WB: 1/1000-1/5000

ICC/IF: 1/100

FC: 1/20

Validations



Western blot detection of Cdk4 in HeLa cell lysates using Cdk4 Rabbit pAb(1:1000 diluted).Predicted band size:34kDa.Observed band size:34kDa.