

CENPH antibody



Catalog Number: 109253

Product name

CENPH antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

CENPH Rabbit Polyclonal antibody. Positive WB detected in NIH/3T3 cells, HeLa cells. Observed molecular weight by Western-blot: 30 kDa and 32 kDa

Preparation

This antibody was obtained by immunization of CENPH recombinant protein (Accession Number: NM_022909). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:500-1:5000

Validations



NIH/3T3 cells were subjected to SDS PAGE followed by western blot with Catalog No:109253(CENPH antibody) at dilution of 1:500