

# MEN1 antibody

Catalog Number: 112653

## Product name

MEN1 antibody

## Specificity

Human; other species not tested.

## Antibody description

MEN1 Rabbit Polyclonal antibody. Positive IP detected in HepG2 cells. Positive WB detected in HepG2 cells, HeLa cells, Jurkat cells. Observed molecular weight by Western-blot: 63-70 kDa

## Preparation

This antibody was obtained by immunization of MEN1 recombinant protein (Accession Number: NM\_001370262). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IP


## Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IP: 1:200-1:2000

## Validations



HepG2 cells were subjected to SDS PAGE followed by western blot with Catalog No:112653(MEN1 antibody) at dilution of 1:500



IP Result of anti-MEN1 (IP:Catalog No:112653, 4ug; Detection:Catalog No:112653 1:500) with HepG2 cells lysate 3200ug.