

# Anti-LRRN3 antibody



Catalog Number: 176653

## Product name

Anti-LRRN3 antibody

## Specificity

Human, Mouse, Rat

## Antibody description

Mouse monoclonal antibody to LRRN3

## Preparation

This antigen of this antibody was recombinant protein within human lrrn3 aa 1-260.

## Formulation

Liquid, 1\*TBS (pH7.4), 0.2% BSA, 50% Glycerol.  
Preservative: 0.05% Sodium Azide.

## Storage

Store at +4°C after thawing. Aliquot store at -20°C. Avoid repeated freeze / thaw cycles.

## Clonality

Monoclonal

## Ig Type

Mouse IgM

## Applications

WB, ICC

## Dilutions

WB: 1:500-1:1,000

ICC: 1:50-1:100

## Validations





Fig1:; Western blot analysis of LRRN3 on different lysates. Proteins were transferred to a PVDF membrane and blocked with 5% BSA in PBS for 1 hour at room temperature. The primary antibody (1/500) was used in 5% BSA at room temperature for 2 hours. Goat Anti-Mouse IgG - HRP Secondary Antibody (HA1006) at 1:5,000 dilution was used for 1 hour at room temperature.; Positive control:; Lane 1: SHSY5Ycell lysate; Lane 2: Rat testis



tissue lysate

Fig2:; ICC staining of LRRN3 in F9 cells (green). Formalin fixed cells were permeabilized with 0.1% Triton X-100 in TBS for 10 minutes at room temperature and blocked with 1% Blocker BSA for 15 minutes at room temperature. Cells were probed with the primary antibody (1/100) for 1 hour at room temperature, washed with PBS. Alexa Fluor®488 Goat anti-Mouse IgG was used as the secondary antibody at 1/1,000 dilution. The nuclear counter stain is DAPI (blue).