

UBC9 antibody



Catalog Number: 107653

Product name

UBC9 antibody

Specificity

Human, Mouse; other species not tested.

Antibody description

UBC9 Mouse Monoclonal antibody. Positive WB detected in 293 cells. Observed molecular weight by Western-blot: 18+54

Preparation

This antibody was obtained by immunization of UBC9 recombinant protein (Accession Number: BC004437). Purification method: Protein A purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Monoclonal

Ig Type

Mouse IgG2a

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:500-1:5000

Validations



HEK-293 cells were subjected to SDS PAGE followed by western blot with Catalog No:107653(UBC9 antibody) at dilution of 1:1000