

CXorf61 antibody



Catalog Number: 109653

Product name

CXorf61 antibody

Clonality

Polyclonal

Immunogen

[Human CXorf61 Recombinant protein \(His tag\)](#)

Ig Type

Rabbit IgG

Specificity

Human, Mouse; other species not tested.

Applications

ELISA, WB

Antibody description

CXorf61 Rabbit Polyclonal antibody. Positive WB detected in mouse brain tissue. Observed molecular weight by Western-blot: 13 kDa

Dilutions

Recommended Dilution:
WB: 1:200-1:2000

Preparation

This antibody was obtained by immunization of CXorf61 recombinant protein (Accession Number: NM_001017978). Purification method: Antigen affinity purified.


Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Validations



mouse brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:109653(CXorf61 Antibody) at dilution of 1:600