

# Anti-MMP12 antibody



Catalog Number: 142353

## Product name

Anti-MMP12 antibody

## Immunogen

[Human MMP-12/MMP12 \(catalytic domain\) recombinant protein](#)

## Specificity

Human MMP-12 / MMP12

## Antibody description

Rabbit Polyclonal to Human MMP12

## Preparation

Produced in rabbits immunized with purified, recombinant Human MMP-12 / MMP12 . MMP-12 / MMP12 specific IgG was purified by Human MMP-12 / MMP12 affinity chromatography.

## Formulation

0.2 µm filtered solution in PBS

## Storage

This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C.

Preservative-Free.

Sodium azide is recommended to avoid contamination (final concentration 0.05%-0.1%). It is toxic to cells and should be disposed of properly. Avoid repeated freeze-thaw cycles.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


IHC-P, ICC/IF, IF

## Dilutions


IHC-P: 0.1-2 µg/mL

ICC/IF: 0.5-3 µg/mL

## Validations



MMP-12 / MMP12 Antibody, Rabbit PAb, Antigen Affinity Purified, Immunohistochemistry



MMP-12 / MMP12 Antibody, Rabbit PAb, Antigen Affinity Purified, Immunofluorescence