

Interferon alpha antibody



Catalog Number: 107353

Product name

Interferon alpha antibody

Specificity

Human; other species not tested.

Antibody description

Interferon alpha Mouse Monoclonal antibody.
Positive WB detected in Recombinant protein.
Observed molecular weight by Western-blot: 22 kDa

Preparation

This antibody was obtained by immunization of Interferon alpha recombinant protein (Accession Number: NM_024013). Purification method: Protein A purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Monoclonal

Ig Type

Mouse IgG2a

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:5000-1:50000

Validations



Recombinant protein were subjected to SDS PAGE followed by western blot with Catalog No:107353(IFNA1 Antibody) at dilution of 1:40000