

Anti-IL18RAP antibody



Catalog Number: 100753

Product name

Anti-IL18RAP antibody

Immunogen

[Rhesus IL18RAP \(His Tag\) recombinant protein](#)

Specificity

Rhesus IL18RAP

Antibody description

Mouse monoclonal to IL18RAP

Preparation

This antibody was produced from a hybridoma resulting from the fusion of a mouse myeloma with B cells obtained from a mouse immunized with purified, recombinant Rhesus IL18RAP (XP_001108108.1; Met1-Arg356). The IgG fraction of the cell culture supernatant was purified by Protein A affinity chromatography.

Formulation

0.2 µm filtered solution in PBS with 5% treha

Storage

This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C.

Preservative-Free.

Sodium azide is recommended to avoid contamination (final concentration 0.05%-0.1%). It is toxic to cells and should be disposed of properly. Avoid repeated freeze-thaw cycles.

Clonality

Monoclonal

Ig Type

Mouse IgG1


Applications

WB

Dilutions

WB: 10-30 µg/ml

Validations



Lanes	A
Items	
Sample (whole cell lysate)	HL-60
Sample Volume (µg/lane)	30
Gel	13% SDS-PAGE reducing gel
Recommended Concentration	10-30 µg/ml
Secondary Antibody	Dylight 800-labeled Antibody to Mouse IgG (H+L) at 1:7500 dilution. Developed using Odyssey imaging system.
Explanation	Predicted band size : 22 kDa Observed band size : 22 kDa

IL18RAP Antibody, Mouse MAb, Western blot